

HAEGARDA® Reconstitution Guide

HAEGARDA® and MIX2VIAL®
(a needleless transfer device) work together to
help you with self-administration

1 CLEAN STOPPERS

Remove the flip caps from both vials (HAEGARDA® and diluent). Wipe rubber stoppers with an antiseptic wipe and allow to dry.



3 PREPARE DILUENT VIAL

Place the diluent vial on an even, clean surface and grip it tightly. Take the Mix2Vial® together with its package and push the spike of the blue adapter end straight down through the diluent vial stopper.



5 TRANSFER DILUENT INTO HAEGARDA® VIAL

Place the HAEGARDA® vial on a firm, even surface. Invert the diluent vial with the Mix2Vial® set attached, and push the spike of the transparent end of the Mix2Vial® straight through the HAEGARDA® vial stopper. The diluent will automatically flow into the HAEGARDA® vial.



2 OPEN THE MIX2VIAL®

Open the Mix2Vial® package by peeling off the lid, leaving the Mix2Vial® inside the package to maintain sterility. Do not remove the Mix2Vial® from the blister package!



4 REMOVE THE MIX2VIAL® PACKAGING

While holding the diluent vial, carefully remove the package from the Mix2Vial® set by pulling it vertically upwards. Make sure to pull away only the package and not the Mix2Vial® set.



BEFORE YOU BEGIN, ASSEMBLE THE SUPPLIES PROVIDED IN THE PRODUCT PACKAGE ON A FLAT SURFACE (INFUSION MAT):

- HAEGARDA® and diluent vials (must be at room temperature)
- Mix2Vial® transfer device
- Sterile syringe
- Subcutaneous infusion set or hypodermic needle
- Alcohol or disinfectant wipes (not included in product package)

Thoroughly clean your working surface using alcohol or disinfectant wipes, and thoroughly wash and dry your hands.

6 DISSOLVE HAEGARDA®

With the diluent and HAEGARDA® vials still attached to the Mix2Vial® set, gently swirl the HAEGARDA® vial to ensure that the powder is fully dissolved. Do not shake the vial.



7 UNSCREW EMPTY DILUENT VIAL

Grasp the transparent HAEGARDA® side of the Mix2Vial® set with one hand, and grasp the blue diluent side with the other hand. Unscrew the set carefully into two pieces. Discard the diluent vial with the blue Mix2Vial® adapter attached.



8 LOAD THE SYRINGE

Draw air into an empty, sterile syringe. With the HAEGARDA® vial upright, screw the syringe into the attached Mix2Vial® set. Inject air into the HAEGARDA® vial.



9 DRAW SOLUTION INTO SYRINGE

Keeping the syringe plunger pressed, turn the system upside down and draw the solution into the syringe by pulling the plunger back slowly.



10 PREPARE ADMINISTRATION SET

Once the solution has been transferred into the syringe, firmly hold the barrel of the syringe (keeping the plunger facing down) and unscrew it from the Mix2Vial® set. The reconstituted solution should be colourless, clear and free from visible particulates or is discoloured. As instructed by a healthcare provider, attach a hypodermic needle or subcutaneous infusion set. Prime the needle or tubing as required and instructed.

